

Step-by-step guide to use PTcalc to interpret the 6MWT and the 2MWT

1. Go to play store link:


https://play.google.com/store/apps/details?id=com.mirza.app to download the apps into your mobile phone




Once the app is downloaded into your mobile phone, update your profile with the following information: Name, company, Current position, Email, Nationality and Institution




3. Then click on 2MWD to use the 2MWT calculator



4. You be directed to this page. This page will show all data that you have saved so far. Click create record to use the calculator



Fill in the required information with regard to your patient's 2MWT. These include patient details for identification (ID/MRN), Age, Gender, Resting Heart Rate, Peak Heart Rate and Distance Walked



The apps will now calculate your estimated 2MWD and the % of expected distance. In this instance, for a female aged 45 years, walked for 120m during her 2MWT with resting HR of 90bpm and peak HR of 120bpm, by using formula 1, she walked only 68% of someone her age without disease could have walked. You can then save your patient data in your database.

